

BATHROOM SOLUTIONS

INNOVATIVE BATHROOM PRODUCTS WITH ELEGANT SCANDINAVIAN DESIGN





"At Ropox we believe it to produce great designs

C

makes sense that can help ther people."

opox offers flexible bathroom solutions to compensate for reduced physical ability. Products are ergonomical, strong and provide extra support during washing, toiletting and bathing.

We are proud of our goals:

- To help people with reduced physical ability to be more independent and to improve, everyday life and self-esteem.

- To make investment in the product economically attractive – in terms of both space-saving (fewer square metres) and fewer care resources.
- To give carers a better working day with a better working environment and fewer days off sick.

That almost completes the story, but we have added that finishing touch:

"At Ropox we believe it makes sense to produce great design, that can help other people". At Harrit & Sørensen we spend a long time in the field to better understand our users. We look for precisely the right angle in a product – we dig for gold. Our philosophy is that design must make a difference; it has to be necessary, arouse emotions and must last - sometimes an entire lifetime."

Nikolaj Sørensen – Partner at Harrit & Sørensen

- 4 ROPOX BATHROOM SOLUTIONS
- 6 ROPOX SWING WASHBASIN & DOCK-IN
- 8 ROPOX GRAB RAILS
- 9 ROPOX SUPPORT WASHBASIN
- 10 ROPOX TOILET SUPPORT ARMS
- 12 ROPOX SHOWER SEAT WITH ARM- AND BACK SUPPORT
- 14 ROPOX CHANGING & SHOWER BED



ROPOX BATHROOM SOLUTIONS



The washbasin with built-in handles, provides support and increased security in the bathroom.



The wash basin gives improved hygiene, as washing of hands can take place before a transfer is made.



Complete bathroom solution with swing washbasin

In bathrooms the swing washbasin, a wheelchair user can easily use the sink sitting down and thereby minimize any risk of falling.

The swing washbasin can be moved as required, allowing more room to maneuver a wheelchair.


The flexible washbasin provides more space to design complete bathrooms with fewer m².


Complete bathroom solution with support washbasin

The support washbasin is robust and particularly suitable for disabled.

Along with the shower seat and the custom-designed support arms, the user has support and the opportunity to sit safely while bathing or using the bathroom in general.

The washbasin can be height adjusted to provide an exact fit for the individual user.





Interior with Swing washbasin

Interior with Support washbasin


ROPOX SWING WASHBASIN & DOCK-IN

The shape and size of the washbasin ensures that seated and standing users can use the basin from the front or sides. A model of the washbasin can also be height adjusted to suit individual users.

- 1. Increased independence Users will be able to use the washbasin and the toilet on their own, giving them greater independence while carers can deal with other tasks.
- **2. Space optimisation** The basin's ability to swing 180° gives more space for wheelchair users and can make smaller rooms more accessible.
- **3. Improved hygiene** Users are able to wash their hands after toiletting before transfer.



VARIANTS



ROPOX SWING WASHBASIN

MATERIALS

PUR FOAM ALUMINIUM STAINLESS STEEL SURFACE TREATED STEEL

FUNCTIONAL DIMENSIONS

THE ARMS SWING 180° THE WASH BASIN SWINGS 180° HEIGHT ADJUSTMENT RANGE 15 CM


ROPOX DOCK-IN

MATERIALS SURFACE TREATED STEEL STAINLESS STEEL COMPACT LAMINATE

LOAD BEARING CAPACITY SWING WASHBASIN & DOCK-IN MAX. 300 KG (150 KG WITH REF. TO DS/EN 12182)

FUNCTIONAL DIMENSIONS HEIGHT ADJUSTMENT RANGE 15 CM

ACCESSORY KIT (INCL. WALLMOUNT/WITHOUT WALLMONT) #40-41150/#40-41145



PLUG WITH OVERFLOW #40-40935




- **1. Useful space** Use of a dock-in means that there is space to put things down within reach.
- Great security for user The use of Dock-in provides a sense of security for users who need extra support.

ROPOX GRAB AND SHOWER RAILS

The handles are a great help for people who need a little more support in their bathroom. They can also serve as towel rails.

- 1. **Greater safety** Provide peace of mind for users and can make self-reliant users independent of carers.
- 2. Easy to clean The smooth surfaces making it easy to clean.





ROPOX SUPPORT WASHBASIN

MATERIALS GELCOATED KOMPOSIT STAINLESS STEEL PLASTIC

LOAD BEARING CAPACITY MAX 400 KG (200 KG WITH REF. TO DS/EN 12182)

FUNCTIONAL DIMENSIONS HEIGHTADJUSTABLE 300 MM.


- BASIC #40-44012 / 40-44013



- FLEXI & ELECTRIC #40-44014 / 40-44015 - #40-44016 / 40-44017



- FIXED #40-44010 / 40-44011



ACCESSORY

CAPE #40-44021 / 40-44022





Basin that is height-adjustable and ideal for bathrooms used by multiple users with individual needs. It may be in hospitals and nursing homes, but also in private homes.

1. Easy to clean - A perfect and elegant design that minimizes dirt from accumulating.


2. Manual/Electric height adjustment -Easy to use without tools, but just a Ropox crank or pushbutton.

ROPOX TOILET SUPPORT ARMS



VARIANTS


10







Random positioning-TheHeight-adjustable toilet support arms arms may be folded up to The height of all toilet support arms can any angle making them more be adjusted by using the wall fitting. convenient to pull down.



Two built-in height settings – With 90 mm. differential, the toilet support arms are extremely useful when standing up and transferring.



Good grip – Their oval shape provides a good, comfortable grip and the colour also makes it ideal for the visually impaired.



Standing/sitting attachment - Makes for a good grip when transferring.



ACCESSORY

STANDING / SITTING ATTACHMENT





TOILET **ROLL HOLDER** #40-40910



MATERIALS

#40-40906

ALUMINIUM PLASTIC STAINLESS STEEL FITTINGS

LOAD BEARING CAPACITY MAX 400 KG (WITH REF. TO DS/ EN 12182, PLS. SEE PICTURES ON PAGE 10)



User-friendly adjustment - The toilet roll holder can be clicked on and rotated 90 degrees to the best position for the user.









ROPOX SHOWER SEAT WITH ARM- & BACK SUPPORT





WITH LEG #40-43011












The Ropox shower seat is ideal for bathrooms where support is needed during bathing and where space is restricted. The shape ensures secure and pleasant showering.


- 1. Heavy load The standard max. load is 150 kg, and the seat is also available with leg to increase the load to 200 kg.
- 2. **Space-saving** The seat can be folded up to give more free space.
- **3. Greater comfort** The material used is hard-wearing and insulating, and the seat does not feel cold or hard.
- **4. Hygiene in focus** The seat is totally sealed polyeurethane which feels warm and has soft round corners and soft edges to facilitate cleaning.



Variable seat height – The variant with leg can be mounted for sitting heights between 420 mm to 540 mm using the discrete adapter




Extremely comfortable – The seat is designed in a comfortable material and in a size and shape that provides great comfort. It also has an easy-tograsp edge to help transfers.




MATERIALS


ALUMINIUM PLASTIC STAINLESS STEEK FITTINGS

LOAD BEARING CAPACITY MAX 500 KG (200 KG WITH A LEG AND 150 KG WITHOUT A LEG WITH REF. TO DS/EN 12182)









ROPOX CHANGING AND SHOWERING BED




Comfortable backrest - giving good support and comfort





Electrical height adju-– Saves time when transferring and makes for a better work environment for carers.

stable bed 300-1000 mm soft and a hard side which can be used as needed. The bed's easy-to-clean design features removable slats.



#40-25038




Added safety - End guards and side rails can be fitted to provide extra safety for users. The see-through bed guard can be folded underneath the shower bed when not needed and for transfers.

MATERIALS PUR SLATS VACUUM FORMED PLASTIC ALUMINIUM SURFACE TREATED STEEL STAINLESS STEEL

LOAD CAPACITY MAX 150 KG WITH REF. TO DS/EN 12182

FUNCTIONAL DIMENSIONS HEIGHT ADJUSTMENT 300-1000 MM = 700 MM

#40-25023




ACCESSORY

#40-25031

In situations where a person needs to be bathed lying down it is important that the bed can be lowered for the person to be transferred - e.g. to a wheelchair, as well as raised to a suitable working height for the carers. Optimum comfort and safety are also important.

Space-saving – The bed can be folded up for additional space.


#40-25036



ABOUT ROPOX

For more than 50 years, Ropox has developed a wide range of products to help people. Ropox develops flexible assistive devices of simple design which look great while improving the daily activities of the users and their helpers.

For further information, please contact Ropox or one of our distributors – www.ropox.com.



ROPOX BATHROOM RANGE

INNOVATIVE AND EXCITING PRODUCTS WHICH HARMONIZE BEAUTIFULLY

ROPOX KITCHEN RANGE

FLEXIBLE AND HEIGHT-ADJUSTABLE KITCHEN SOLUTIONS

ROPOX TABLE RANGE

HEIGHT-ADJUSTABLE TABLES ADAPTED TO INDIVIDUAL NEEDS Copyrights © All Rights Reserved Ropox A/S 2014

taylordolman

tel: 0845 838 0600 email: info@taylordolman.com www.taylordolman.com

