"An electrical shower/changing bed with a high degree of comfort"





OPTIMUM WORKING CONDITIONS FOR THE CARE WORKERS AND A SAFE AND COMFORTABLE EXPERIENCE FOR THE USERS

Lying on comforable slats

The shower bed's design ensures that the user always lies comfortably on the slats and that the water will drain straight away so that the user is not in contact with their dirty water – this is collected in the drip tray below the lying surface of the slats and taken away by flexible waste pipe.

The slats can be lifted out individually for easy cleaning between users. The reversible slats have a soft, warm and pleasant feel. They bend slightly downwards on one side, so the user lies safely and securely in the shower, and water and soap can run quickly into the tray and away. The other side bends slightly upwards, which makes it easier for the care workers to turn a user when helping them change.

Security in focus


The height is adjusted with a hand control, providing the care workers with a healthy work posture and easy access to the shower bed. It also makes it more manageable to help a user onto the shower bed from, e.g., a wheelchair. The mobile changing/shower bed can be supplied with a battery solution, so you can move it without considering if cords or cables get in the way. It is also possible to mount a central brake system so the carers can brake the shower bed easily and quickly. All in all, the Mobilio shower/changing bed is a flexible, comfortable and safe solution for both users and care workers.

Easy transportation

The height adjustment gives the carers the right working height, and it reduces the amount of person lifting and moving because it is easy to help a user from, e.g., the bedroom to the bathroom. The users will find that changing and showering becomes a safe and comfortable experience.



ROPOX MOBILIO SHOWER BED



- * The shower bed is easy to adjust in height with the electrical hand control.
- * The lying area is 18 cm shorter than the length of the shower bed.
- * It is possible to mount end guard (optional) at both ends.
- * Bed guards for the sides are included.
- * If fitted with a battery (optional), the product can be used also when charged.

ACCESSORIES

Back support Item no.: 40-25031



Head rest Item no.: 40-25090



End guard Item no.: 40-25038



Central brake Length 146 cm - Item no.: 40-30753 Length 178 cm - Item no.: 40-30760 Length 210 cm - Item no.: 40-30763



Battery Cover for Item no.: 40-30762 Length 146



Cover for bed guard Length 146 cm ltem no.: 40-25782 Length 178 cm, 210 cm ltem no.: 40-25783



Find more information about measurements, technical data, installation interfaces and much more under the product and in the instruction manuals on our website.



www.taylordolman.com

