

Etac Molift Partner 255

The flagship of Etac Molift's mobile lifters.



*Find optional suspensions
and detailed information
at page 26–27.*

Vertical lifting pattern and excellent range

Etac Molift Partner 255 provides an ideal solution for lifting with high complexity situations in acute and critical care areas in both hospitals and care facilities. With an excellent lifting range: 37 – 179 cm (14.5" – 70.4") the Molift Partner 255 allows easy lifting from the floor as well as high surfaces.

The straight vertical column in combination with the fixed suspension provides excellent stability and maintains the user at a good distance from the column.

Stability, maneuverability and lifting capacity

With a total weight of 42.5 kg (94 lbs) Molift Partner 255 is one of the lightest in its category, and combined with excellent construction design and wheel solution, it becomes very easy to maneuver with a safe working load of 255 kg (562 lbs). Equipped with environmentally friendly NiMH-battery.

Description

Etac Molift Partner 255 excl. suspension
Standard 13000
Low base 13050
Supporting arms 13100

Material

Aluminium, Steel

Service Software

Included

Total weight

Regular: 42.5 kg (94 lbs)
Low base: 45 kg (99 lbs)
With support arms: 50 kg (100.5 lbs)

Height

190 cm (74.8")

Length

122 cm (48.0")

Width chassis/legs outer

70 cm (27.6")

Width open legs inner

107 cm (42.1")

Lifting height

37 – 179 cm
(14.5" – 70.4")

Lifting interval

142 cm (55.9")

Leg height

11.5 cm (4.5")
Low base: 6 cm (2.4")
With support arms: 11.5 cm (4.6")

Battery

NiMH 26.4 V – 2.2 Ah

Battery charge time

Ca 3 hours

Lifting speed

3 cm/ sec
(1.6" / sec)

Rated performance

50 lifts at: 75 kg and 50 cm
(165 lbs and 19.7")

Cleaning instructions

Use pH-neutral cleaning agents. If disinfection is required, use isopropyl alcohol. Do not use abrasive or corrosive substances for cleaning.



Features



High safe working load

Molift Partner 255 has a safe working load of 255 kg (562 lbs)



Low base

Molift Partner 255 is also available with an extra low base chassis. A base height of only 6 cm (2.4") is one of the lowest in the market.



Walking arms

Molift Partner 255 is available with additional Support Arms. This in combination with the maximum height provides an ideal solution for walking training.



4-point suspension

Ensures a comfortable and spacious lifting position. Combined with the vertical lifting pattern it creates a safe and stable transfer.

Accessories



Weight Scale Arm – Molift Partner 255 1331002

The Molift Partner 255 scale is placed on top of the lifting arm, thereby subtracting nothing from the lifting height. The scale is approved Class III and can also calculate BMI (Body Mass Index).



Etac Molift Partner 255 has a vertical lifting pattern which keeps the user at an equidistant length from the column during the entire lift, creating a safe and stable position.