




New, improved lifting hanger from Guldmann

Guldmann is one of the leading and most innovative companies in the market of stationary and mobile lifters. Guldmann's newly developed lifting hanger makes the everyday tasks easier for both user and helper.

The new hanger replaces our previous hangers with 2 and 4 attachments in our product range. Function, size and attachment to the hoist remains the same in order for the hangers to be able to be used in the same systems as today.

The new design gives the hanger multiple benefits:

-  **Lower weight**
The new hanger is 30 % lighter than the previous 4 point hanger.
-  **Increased strength**
The new hanger's maximum load has been increased to 375 kg. Making it usable for the heavier users.
-  **Improved comfort**
The hanger has soft, seamless surfaces that make it easier to clean and more comfortable to touch. Furthermore, the ends of the hanger are equipped with "shock absorbers" in plastic.

	Discontinued item numbers		NEW item numbers	
	2 attachment Points Hanger	4 attachment Points Hanger	4 attachment Points Hanger	4 attachment Points Hanger Adjustable
XS	550058	554260	556870	556910
S	550057	554259	556880	556920
M	550316	554258	556890	556930

| Time to care |



GH lifting hanger (2016) Technical specifications

Functions

Lifting capacity, max. : 375 kg
Suspension points : 4

Weight and Materials

Medium, totally : 1,8 kg
SMALL, total : 1,6 kg
X-SMALL, total : 1,35 kg

Lifting hanger : Rustproof steel covered with polyurethane foam
Suspension points : Rustproof steel

Measurements

556890 MEDIUM

A : 600 mm
B : 534 mm
C : 128 mm

556880 SMALL

A : 500 mm
B : 440 mm
C : 100 mm

556870 X-SMALL

A : 395 mm
B : 343 mm
C : 52 mm

Accessories/Miscellaneous



Lifting sling : See brochure

Labelling

The product is manufactured in compliance with the Council Directive 93/42/EEC of June 14th 1993, including amendments, as medical device class 1.

Classification

according to ISO 9999 : 12 36 03



| Time to care |